

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient Comfort Care Cuffs – Single Hose

Product #	Description	Packaging	Use with:*
M1571A	Comfort Care Cuff, Infant Bladder width: 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm	1 cuff	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) 78352A, 78352C, 78354A, 78354C, M1008A, M1008B, M1176A, M1350B, M1350C, 860335, 860338, 861290, 861493, 861494, 861495, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863275, 863276, 863277, 863278, 863279, 863283, M2702A, M2703A, M2704A, M2705A, M3001A, M3002A, M3535A, M3536A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536M7, M3536M8, M3536M9, M3536MC, M8102A, M8105A, M8105AS, M8105AT
M1572A	Comfort Care Cuff, Pediatric Bladder width: 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm		
M1573A	Comfort Care Cuff, Small Adult Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28cm		
M1573XL	Comfort Care Cuff, Small Adult XL Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28cm		
M1574A	Comfort Care Cuff, Adult Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient Comfort Care Cuffs – Single Hose (continued)

Product #	Description	Packaging	Use with:*
M1574XL	Comfort Care Cuff, Adult XL Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm	1 cuff	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) 78352A, 78352C, 78354A, 78354C, M1008A, M1008B, M1176A, M1350B, M1350C, 860335, 860338, 861290, 861493, 861494, 861495, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863275, 863276, 863277, 863278, 863279, 863283, M2702A, M2703A, M2704A, M2705A, M3001A, M3002A, M3535A, M3536A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536M7, M3536M8, M3536M9, M3536MC, M8102A, M8105A, M8105AS, M8105AT
M1575A	Comfort Care Cuff, Adult Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34cm to 43cm		
M1575XL	Comfort Care Cuff, Large Adult XL Bladder width: 6.3" (16 cm) Bladder length: 13.8" (35 cm) Limb circumference: 34cm to 43cm		
M1576A	Comfort Care Cuff, Large Adult XL Bladder width: 7.9" (20 cm) Bladder length: 16.9" (43 cm) Limb circumference: 42cm to 54cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient Easy Care Cuffs - Single Hose

Product #	Description	Packaging	Use with:*
M4552B / M4552B5	Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm	1 cuff / 5 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478
M4553B / M4553B5	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm		
M4554B / M4554B5	Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference: 20.5cm to 28.5cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient Easy Care Cuffs - Single Hose (continued)

Product #	Description	Packaging	Use with:*
M4555B / M4555B5	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference: 27.5cm to 36cm	1 cuff / 5 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
M4556B / M4556B5	Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4557B / M4557B5	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference: 35cm to 45cm		
M4558B / M4558B5	Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4559B / M4559B5	Thigh, grey Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference: 44cm to 56cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient Easy Care Cuffs - Double Hose

Product #	Description	Packaging	Use with:*
M4562B	Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm	1 cuff	For use with double hose monitors, including Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
M4563B	Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm		
M4564B	Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference: 20.5cm to 28.5cm		
M4565B	Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4566B	Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4567B	Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4568B	Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4569B	Thigh, grey. Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference: 44cm to 56cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient NBP Cuff Kits – Single Hose

Product #	Description	Packaging	Use with:*
864288	Easy Care Adult Kit – 4 sizes Small Adult, Adult, Large Adult, Thigh	consists of 4 cuffs: M4554B, M4555B, M4557B, M4559B	Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
864289	Easy Care Pediatric Kit – 4 sizes Infant, Pediatric, Small Adult, Adult	consists of 4 cuffs: M4552B, M4553B, M4554B, M4555B	
864290	Easy Care Assortment Kit – 6 sizes Infant, Pediatric, Small Adult, Adult, Large Adult, Thigh	consists of 6 cuffs: M4552B, M4553B, M4554B, M4555B, M4557B, M4559B	
864291	Easy Care Adult Long Kit – 6 sizes Small Adult, Adult, Adult Long, Large Adult, Large Adult Long, Thigh	consists of 6 cuffs: M4554B, M4555B, M4556B, M4557B, M4558B, M4559B	
M1577A	Comfort Care Pediatric Kit – 4 sizes Adult, Infant, Pediatric, Small Adult	consists of 4 cuffs: M1571A, M1572A, M1573A, M1574A	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Reusable NBP Cuffs



Non-Invasive Blood Pressure – Reusable NBP Cuffs

Multi Patient NBP Cuff Kits – Single Hose (continued)

Product #	Description	Packaging	Use with:*
M1578A	Comfort Care Adult Kit – 4 sizes Small Adult, Adult, Large Adult, Thigh	consists of 4 cuffs: M1573A, M1574A, M1575A, M1576A	78352A, 78352C, 78354A, 78354C, M1008A, M1008B, M1176A, M1350B, M1350C, 860335, 860338, 861290, 861493, 861494, 861495, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863275, 863276, 863277, 863278, 863279, 863283, M2702A, M2703A, M2704A, M2705A, M3001A, M3002A, M3535A, M3536A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536M7, M3536M8, M3536M9, M3536MC, M8102A, M8105A, M8105AS, M8105AT
M1579A	Comfort Care Assortment Kit – 6 sizes Infant, Pediatric, Small Adult, Adult, Large Adult, Thigh	consists of 6 cuffs: M1571A, M1572A, M1573A, M1574A, M1575A, M1576A	
M1579XL	Comfort Care Adult Long Kit-3 sizes Small Adult XL, Adult XL, Large Adult XL	consists of 3 cuffs: M1573XL, M1574XL, M1575XL	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Disposable NBP Cuffs



Non-Invasive Blood Pressure – Disposable NBP Cuffs

Multi Patient Disposable Cuffs - Multi Care

Product #	Description	Packaging	Use with:*
989803183311	Infant (orange) Bladder width = 2.2" (5.5cm) Bladder length = 6.1" (15.4cm) Limb circumference: 10cm to 15cm	10 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478
989803183321	Pediatric (green) Bladder width = 3.1" (8cm) Bladder length = 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm		
989803183331	Small Adult (light blue) Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28cm		
989803183341	Adult (navy blue) Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm		
989803183351	Adult XL (navy blue with edge striping) Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm		
989803183361	Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34cm to 43cm		
989803183371	Thigh (grey) Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42cm to 54cm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Disposable NBP Cuffs (single-patient use)



Non-Invasive Blood Pressure – Disposable NBP Cuffs (single-patient use)

Gentle Care Disposable Cuffs - Double Hose

Product #	Description	Packaging	Use with:*
① M1866S CE	Neonatal, soft, Single-patient cuff, size #1, single hose	40 cuffs per box	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) NBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
① M1868S CE	Neonatal soft Single-patient cuff, size #2, single hose		
① M1870S CE	Neonatal soft Single-patient cuff, size #3, single hose		
① M1872S CE	Neonatal soft Single-patient cuff, size #4, single hose		
① M1873S CE	Infant, soft, Single-patient cuff, size #5, single hose		
① 989803167551 CE	Neonatal, soft, Single-patient, cuffs, single hose consists of: 5 x M1866S, 10 x M1868S, 20 x M1870S, 10 x M1872S & 5 x M1873S)	50 cuffs per box	NBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)



Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)

Single-patient, vinyl, neonatal/infant cuffs

Product #	Description	Packaging	Use with:*
1 M1866B CE	Neonatal Single-patient cuff, size 1 Circumference: 3.1cm to 5.7cm (1.2" to 2.2"), 1 hose. New connector configuration, replaces the M1866A.	40 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) NBP module: M1008B; Multi-Measurement servers: M3000A, M3001A, M3002A; Patient monitors M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082; Legacy HP monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
1 M1868B CE	Neonatal Single-patient cuff, size 2 Circumference: 4.3cm to 8.0cm (1.7" to 3.1"), 1 hose. New connector configuration, replaces the M1868A.		
1 M1870B CE	Neonatal Single-patient cuff, size 3 Circumference: 5.8cm to 10.9cm (2.3" to 4.3"), 1 hose. New connector configuration, replaces the M1870A.		
1 M1872B CE	Neonatal Single-patient cuff, size 4 Circumference: 7.1cm to 13.1 cm (2.8" to 5.2"), 1 hose. New connector configuration, replaces the M1872A.		
1 M1873B CE	Infant, single patient cuff, size 5 Circumference: 10.0cm to 15.0cm, (3.9" to 5.9"), 1 hose New connector configuration.		
1 989803167541 CE	Neonatal Single-patient NBP kit contains 50 cuffs. Incorporates new connector configuration and replaces M1820-60020.	5x M1866B size 1, 10x M1868B size 2, 20x M1870B size 3, 10x M1872B size 4, 5x M1873B size 5	NBP module: M1008B; Multi-measurement servers: M3000A, M3001A, M3002A, Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors; 78352A, 78352C, 78834C, M8105AS

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)



Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)

Single Care Disposable Cuffs

Product #	Description	Packaging	Use with:*
989803182281	Pediatric (green) Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14–21.5cm	20 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478
989803182291	Small Adult (light blue) Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference: 20.5cm to 28.5cm		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
989803182301	Adult (navy blue) Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference: 27.5cm to 36cm		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
989803182311	Adult XL (navy blue with edge striping) Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
989803182321	Large Adult (burgundy) Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference: 35cm to 45cm		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)

Gentle Care Disposable Cuffs – Double Hose

Product #	Description	Packaging	Use with:*
M4582B	Infant (orange) Bladder width 2.2" (5.5cm), Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm	10 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
M4583B	Pediatric (green) Bladder width 3.1" (8cm), Bladder length: 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm		
M4584B	Small Adult (light blue) Bladder width: 4.2" (10.6cm), Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28.5cm		
M4585B	Adult (navy blue) Bladder width: 5.3" (13.5cm), Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4586B	Adult Extra-Long (navy blue) Bladder width: 5.3" (13.5cm), Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		For use with double hose monitors
M4587B	Large Adult (burgundy) Bladder width: 6.7" (17cm), Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4588B	Large Adult Extra-Long (burgundy) Bladder width: 6.7" (17cm), Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4589B	Thigh (grey) Bladder width: 8.3" (21cm), Bladder length: 18.9" (48cm) Limb circumference: 44cm to 56cm		

Non-Invasive Blood Pressure – Disposable NBP Cuffs (single patient use)

Gentle Care Disposable Cuffs – Single Hose

Product #	Description	Packaging	Use with:*
M4572B	Infant (orange) Bladder width 2.2" (5.5cm), Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm	10 cuffs	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Xper Flex Cardio Physiomonitors System: 860335, 860338
M4573B	Pediatric (green) Bladder width 3.1" (8cm), Bladder length: 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm		
M4574B	Small Adult (light blue) Bladder width: 4.2" (10.6cm), Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28.5cm		
M4575B	Adult (navy blue) Bladder width: 5.3" (13.5cm), Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4576B	Adult Extra-Long (navy blue) Bladder width: 5.3" (13.5cm), Bladder length: 11.8" (30cm) Limb circumference: 27.5cm to 36cm		
M4577B	Large Adult (burgundy) Bladder width: 6.7" (17cm), Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4578B	Large Adult Extra-Long (burgundy) Bladder width: 6.7" (17cm), Bladder length: 15" (38cm) Limb circumference: 35cm to 45cm		
M4579B	Thigh (grey) Bladder width: 8.3" (21cm), Bladder length: 18.9" (48cm) Limb circumference: 44cm to 56cm		
			Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Xper Flex Cardio Physiomonitors System: 860335, 860338; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure - NBP Accessories and Hoses



Non-Invasive Blood Pressure – NBP Accessories and Hoses

Blood Pressure Air Hoses

Product #	Description	Packaging	Use with:*
1 989803136931	Adult/Pediatric NBP air hose Length: 11.71ft (3.57m)	1 air hose	VS1 monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061
M1598B	Adult NBP air hose Length: 4.92ft (1.5m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
2 M1599B	Adult NBP air hose length: 10ft (3.0m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; IntelliVue Monitors: MP2: M8102A, MP5: M8105A, M8105AS, M8105AT, 866044, 866045, 866046; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863275, 863276, 863277, 863278, 863279, 863283, 863285, 863286, 863287, 863288, 863289, 862474, 862478; Xper Flex Cardio Physiomonitors: 860335, 860338; Fetal Monitors: M2702A, M2703A, M2704A, M2705A, M1351A, M1353A, M1350B, M1350C
3 M3918A	Adult/Pediatric air hose Length: 10ft (3.0m)		A1: M3921A, M3922A, M3923A, M3924A, A3: M3926A, M3927A, M3928A, M3929A; ICG: 862146; C1 Monitors: 863051, 863052, 863053, 863054
M1596C	Neonatal NBP air hose Length: 4,92ft (1.5m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) M1008B, 863063, 863069, 863275, 863073, M8102A, M8105A, M3001A, M3000A, M3002A
M1597C	Neonatal NBP air hose Length: 10ft (3.0m)	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) M1008B, 863063, 863069, 863275, 863073, M8102A, M8105A, M3001A, M3000A, M3002A	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure - NBP Accessories and Hoses



1



2

Non-Invasive Blood Pressure – NBP Accessories and Hoses

Reusable Value Care Cuffs

Product #	Description	Limb Circumference	Bladder width	Use with:*
989803160861	Large adult	34 cm to 43 cm	16 cm	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302),
989803160851	Adult extra-long	27 cm to 35 cm	13 cm	
989803160841	Adult	27 cm to 35 cm	13 cm	
989803160831	Small adult	20 cm to 28 cm	10.5 cm	
989803160821	Pediatric	14 cm to 21 cm	8 cm	
989803160811	Infant	10 cm to 15 cm	5.5 cm	

NBP Cuff Caps

Product #	Description	Packaging	Use with:*
① 989803151601	Cuff caps for single tube cuffs	10 caps	n/a
② 989803151611	Cuff caps for double tube cuffs	20 caps	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure - NBP Accessories and Hoses



1



2

Non-Invasive Blood Pressure – NBP Accessories and Hoses

Bulb and Valve Kit for Manual Sphygmomanometer

Product #	Description	Packaging	Use with:*
1 M4580A	Bulb and Valve Kit Manual Sphygmomanometer Replacement coil tubing (M4581A) is available for this product	1 bulb and valve	2-hose reusable cuff models M456xB and soft single patient disposable models M458xB only.
2 M4581A	Coiled Tubing Bulb and Valve Kit, Replacement, length: 8"	1 tubing	M4580A Bulb and Valve Kit

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Non-Invasive Blood Pressure - NBP Accessories and Hoses



1



2



3



4



5



6



7



8



9



10



11



12

Non-Invasive Blood Pressure – NBP Accessories and Hoses

NBP Cuff Adapters

Product #	Description	Packaging	Use with:*
1 M4550A	1 Hose adapter (Welch-Allyn) Threaded female style connector	1 adapter	Philips 1-hose reusable and disposable cuffs excluding 40401x series
2 M4551A	2 Hose Y-adapter (Welch-Allyn) Y style cuff adapter		Philips 2-hose reusable and disposable cuffs
3 M4560A	1 Hose adapter Female quick-connect style connector		Philips 1-hose reusable and disposable cuffs excluding 40401x series
4 M4561A	2 Hose Y-adapter Female quick-connect style connector		Philips 2-hose reusable and disposable cuffs
5 M4570A	1 Hose adapter (Draeger) Luer-lock style cuff adapter		Philips 1-hose reusable and disposable cuffs excluding 40401x series
6 M4571A	2 Hose Y-Adapter (Draeger) Y Style		Philips 2-hose reusable and disposable cuffs
M4590A	1 Hose Adapter (Colin, Datascope-Passport, Spacelabs, Welch-Allyn, Mindray) Locking luer style connection		Philips 1-hose reusable and disposable cuffs excluding 40401x series
7 M4591A	1 Hose Adapter (Welch-Allyn Protocol) Male quick-connect style connector		Philips 2-hose reusable and disposable cuffs
8 M4592A	2 Hose Adapter (GE/Marquette) Male quick connect style connector		
9 M4593A	2 Hose Adapter (Datex/Ohmeda) One male and one female quick-connect style connector		
10 M4594A	Y Adapter (Spacelabs, Siemens, Mindray) Bayonet style connector		
11 M4595A	Y Adapter (Colin, Datascope-Passport, Spacelabs, Welch-Allyn, Mindray) Locking luer style connector		
12 M4596A	Y Adapter (Welch-Allyn Protocol) Male quick-connect style connector		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.